


The Tecres preformed spacer for hip has a load-bearing structure in stainless steel coated with antibiotic bone cement and it resembles a femoral prosthesis.

**Spacer-G:** loaded with Gentamicin

(3 sizes of head with long stem and 3 sizes of head with short stem)

**Spacer-G Flat Stem:** loaded with Gentamicin with flat stem designed to preserve the greater trochanter

(3 sizes of head with long stem and 3 sizes of head with short stem)



**SPACER G - short stem**

With gentamicin

Available in 3 sizes of head: 46/54/60 mm



**SPACER G FLAT STEM - short stem**

With gentamicin

Available in 3 sizes of head: 46/54/60 mm



**SPACER G - long stem**

With gentamicin

Available in 3 sizes of head: 46/54/60 mm



**SPACER G FLAT STEM - long stem**

With gentamicin

Available in 3 sizes of head: 46/54/60 mm